


Tissue

Sample Tags

Integrated Tissue Sampling and Animal ID

We introduced the industry-leading **one**-piece, Z1 No-Snag-Tag, and now we're bringing you the convenience of tissue sampling and animal ID application in **one** simple step with the Tissue Sample Tags. The cross-functional Tissue Tagger seamlessly installs RFID and Tissue Sample Tags quickly and cleanly for a simple **one**-step process. Dry desiccant vial samples offer high-quality preservation and can be stored at room temperature for up to six months. This process minimizes the risk of sample destruction, ensuring integrity for genetic and disease testing.

- Visual Tag guaranteed for the life of the animal
- Quick, clean collection for genetics and disease testing
- Accurate offspring pairing for improved reproduction results
- Increased herd credibility
- Prevents recording errors within the herd
- Matching barcodes with Tissue Sample, RFID and visual tags
- Store sample up to 6 months at room temperature
- Custom laser printing available
- Certified PAS44



Herdsmen Tags



Stockman Tags



Z1 No-Snag-Tags



REQUIRES
Z1 NO-SNAG-TAGGER

Tissue Tagger



800 0800-123

EID



Tag Marking Pen



700 2500-239



◀◀ Learn More

VISIT [ZTAGS.COM](https://www.ztags.com)

Customer Support **800-433-3112** or customerservice@datamars.com

Z Tags is registered trademark of Datamars, Inc. © 2024 Datamars, Inc. – All rights reserved. REV 07182024



TSO SPECIFICATIONS

All TSOs are yellow, laser printed and applied using the Tissue Tagger. Matching numbers available with EID or Visual Tag.

Ea = Each; Pk = Pack

	SKU	Type of Tag	Matched Set	EID	USDA	Quantity
TSO (Box of 10)	700 2001-391	1pc Tissue Sample	—	—		Box of 10 TSO
TSO + ComfortEar Slim F4Z M2Z HDX 840	700 2510-604	1pc Tissue Sample, 2pc EID	Yes	HDX White	840 Official	Box of 10 TSO; 10 EIDs/set
TSO + ComfortEar Slim F4Z M2Z FDX 840	700 2510-606	1pc Tissue Sample, 2pc EID	Yes	FDX Yellow	840 Official	Box of 10 TSO; 10 EIDs/set
TSO + ComfortEar Slim F4Z M2Z HDX 964	700 2510-608	1pc Tissue Sample, 2pc EID	Yes	HDX White	964 Regulated	Box of 10 TSO; 10 EIDs/set
TSO + ComfortEar Slim F4Z M2Z FDX 964	700 2510-610	1pc Tissue Sample, 2pc EID	Yes	FDX Yellow	964 Regulated	Box of 10 TSO; 10 EIDs/set

TST SPECIFICATIONS

All TSTs are yellow, custom laser printed and applied using the Tissue Tagger. Matching numbers available with EID or Visual Tag.

Ea = Each; Pk = Pack

	SKU	Type of Tag	Matched Set	EID	USDA	Quantity
TST + Laser One Side	700 2501-589	2pc Tissue Sample VID	—	—		1 Pk {4ea}
TST + Laser Two Sides	700 2501-590	2pc Tissue Sample VID	—	—		1 Pk {4ea}
TST + ComfortEar Slim F4Z HDX 840; Laser 1 Side	700 2510-612	2pc Tissue Sample VID; 2pc EID	Yes	HDX White	840 Official	1 Pk {4ea}; 4 EIDs/set
TST + ComfortEar Slim F4Z FDX 840; Laser 1 Side	700 2510-614	2pc Tissue Sample VID; 2pc EID	Yes	FDX Yellow	840 Official	1 Pk {4ea}; 4 EIDs/set
TST + ComfortEar Slim F4Z HDX 964; Laser 1 Side	700 2510-616	2pc Tissue Sample VID; 2pc EID	Yes	HDX White	964 Regulated	1 Pk {4ea}; 4 EIDs/set
TST + ComfortEar Slim F4Z FDX 964; Laser 1 Side	700 2510-618	2pc Tissue Sample VID; 2pc EID	Yes	FDX Yellow	964 Regulated	1 Pk {4ea}; 4 EIDs/set
TST + ComfortEar Slim F4Z HDX 840; Laser 2 Sides	700 2510-620	2pc Tissue Sample VID; 2pc EID	Yes	HDX White	840 Official	1 Pk {4ea}; 4 EIDs/set
TST + ComfortEar Slim F4Z FDX 840; Laser 2 Sides	700 2510-622	2pc Tissue Sample VID; 2pc EID	Yes	FDX Yellow	840 Official	1 Pk {4ea}; 4 EIDs/set
TST + ComfortEar Slim F4Z HDX 964; Laser 2 Sides	700 2510-624	2pc Tissue Sample VID; 2pc EID	Yes	HDX White	964 Regulated	1 Pk {4ea}; 4 EIDs/set
TST + ComfortEar Slim F4Z FDX 964; Laser 2 Sides	700 2510-626	2pc Tissue Sample VID; 2pc EID	Yes	FDX Yellow	964 Regulated	1 Pk {4ea}; 4 EIDs/set